

Product Datasheet

HIST1H3A (phospho-S10) antibody (orb214781)

Description

Rabbit polyclonal antibody to HIST1H3A

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat, Zebrafish

Conjugation

Unconjugated

Tested

IF, IH, WB

Applications

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S10 of human Histone H3 protein. The exact sequence is proprietary.

Target

HIST1H3A; HIST1H3B; HIST1H3C; HIST1H3D; HIST1H3E; HIST1H3F; HIST1H3G; HIST1H3H; HIST1H3I; HIST1H3J; HIST2H3A; HIST2H3C; HIST2H3D; H3F3A; H3F3B

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P68431, P68433, Q6LED0, P84244, P84228, Q71DI3, P84245, P84243

Entrez

8352, 8356, 8354, 8968, 15078, 8357, 319152, 8358, 8355, 291159, 100361558, 8350, 15077, 8353, 8351

Dilution Range


WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500

Expiration Date

12 months from date of receipt.



Western blot analysis of Histone H3 (Pho...



Immunohistochemical analysis of Histone ...



Immunofluorescent analysis of Histone H3...