

## Product Datasheet

### PRDX4 ELISA Kit (Rat) (orb2148256)

|                            |                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------|
| <b>Description</b>         | PRDX4 ELISA Kit (Rat)                                                               |
| <b>Range</b>               | 0.156-10ng/mL                                                                       |
| <b>Tested Applications</b> | ELISA                                                                               |
| <b>Storage</b>             | Store as indicated in product manual.                                               |
| <b>Note</b>                | For research use only                                                               |
| <b>Application notes</b>   | Assay Info: Assay Methodology: Quantitative Sandwich ELISA                          |
| <b>Sample Types</b>        | Serum, Plasma, Tissue homogenates, Cell culture supernates, Other biological fluids |
| <b>MW</b>                  | 31 kDa                                                                              |
| <b>Uniprot ID</b>          | <b>Q9Z0V5</b>                                                                       |
| <b>Sensitivity</b>         | 0.078 ng/mL                                                                         |
| <b>NCBI</b>                | <b>NP_445964.1</b>                                                                  |
| <b>Expiration Date</b>     | Please enquire.                                                                     |