

Product Datasheet

RGS3 ELISA Kit (Mouse) (orb2148304)

Description	RGS3 ELISA Kit (Mouse)
Range	78-5000pg/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue homogenates, Cell culture supernates, Other biological fluids
MW	106 kDa
Uniprot ID	Q9DC04
Sensitivity	39.8 pg/mL
NCBI	NP_001075119.1
Expiration Date	Please enquire.