

Product Datasheet

Hmox1 ELISA Kit (Mouse) (orb2148583)

Description	Hmox1 ELISA Kit (Mouse)
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue homogenates, Cell culture supernates, Other biological fluids
MW	33 kDa
Uniprot ID	P14901
Sensitivity	0.078 ng/mL
NCBI	NP_034572.1
Expiration Date	Please enquire.