

Product Datasheet

TCF3 ELISA Kit (Human) (orb2148933)

Description	TCF3 ELISA Kit (Human)
Range	0.312-20ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Competitive ELISA
Sample Types	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
MW	68 kDa
Uniprot ID	P15923
Sensitivity	0.19 ng/mL
NCBI	NP_001129611.1
Expiration Date	Please enquire.