

Product Datasheet

SPAG9 ELISA Kit (Human) (orb2149119)

Description	SPAG9 ELISA Kit (Human)
Range	0.156-10ng/ml
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
MW	146 kDa
Uniprot ID	O60271
Sensitivity	0.056 ng/mL
NCBI	NP_001123999.1
Expiration Date	Please enquire.