

Product Datasheet

SORBS1 ELISA Kit (Human) (orb2149597)

Description	SORBS1 ELISA Kit (Human)
Range	0.312-20ng/ml
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
MW	11 kDa
Uniprot ID	Q9BX66
Sensitivity	0.12 ng/mL
NCBI	NP_001030126.1
Expiration Date	Please enquire.