

Product Datasheet

HUWE1 ELISA Kit (Human) (orb2149658)

Description	HUWE1 ELISA Kit (Human)
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
MW	481 kDa
Uniprot ID	Q7Z6Z7
Sensitivity	0.067 ng/mL
NCBI	NP_113584.3
Expiration Date	Please enquire.