

Product Datasheet

ATF3 ELISA Kit (Mouse) (orb2151305)

Description	ATF3 ELISA Kit (Mouse)
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue homogenates, cell lysates
MW	21 kDa
Uniprot ID	Q60765
Sensitivity	0.063ng/mL
NCBI	NP_031524.2
Expiration Date	Please enquire.