



Product Datasheet

H2BFS antibody (orb215195)


Description	Rabbit polyclonal antibody to H2BFS
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Zebrafish
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Histone H2B. The exact sequence is proprietary.
Target	H2BFS
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P57053
Entrez	102724334
Dilution Range	WB: 1-500-1-1000, IHC-P: 1-100-1-200, IF/ICC: 1-100-1-500, ChIP: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2B exp...



Immunohistochemical analysis of Histone ...



Immunofluorescent analysis of Histone H2...