

Product Datasheet

HIF1AN ELISA Kit (Human) (orb2152002)

Description	HIF1AN ELISA Kit (Human)
Range	31.25-2000pg/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NWT6
Sensitivity	14.4pg/mL
NCBI	NP_060372.2
Expiration Date	Please enquire.