


Product Datasheet

NBPF12 antibody (orb215219)

Description	Rabbit polyclonal antibody to NBPF12
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human NBPF-pan. The exact sequence is proprietary.
Target	NBPF12
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q3BBV0 , Q5TAG4 , Q8N660 , Q6P3W6 , Q5TI25 , Q3BBV1 , Q5SXJ2 , Q3BBW0
Entrez	284565 , 400818 , 25832 , 55672
Dilution Range	WB: 1-500-1-1000, IP: 1-10-1-100
Expiration Date	12 months from date of receipt.



Western blot analysis of NBPF-pan expres...



Immunohistochemical analysis of NBPF-pan...