



Product Datasheet

Apolipoprotein AI antibody (orb215363)

Description	Rabbit polyclonal antibody to Apolipoprotein AI
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human APOA1BP. The exact sequence is proprietary.
Target	APOA1BP
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8K4Z3 , Q8NCW5 , B0BNM1
Entrez	246703 , 128240 , 295229
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.




Western blot analysis of APOA1BP express...



Immunohistochemical analysis of APOA1BP

...



Immunofluorescent analysis of APOA1BP st...