

Product Datasheet

CAV1 ELISA Kit (Chicken) (orb2156265)

Description	CAV1 ELISA Kit (Chicken)
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	2°C to 8°C, -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P35431
Sensitivity	0.064 ng/mL
NCBI	NP_001099134.1
Expiration Date	Please enquire.