

Product Datasheet

ANGPT2 ELISA Kit (Rabbit) (orb2158205)

Description	ANGPT2 ELISA Kit (Rabbit)
Range	0.312-20ng/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich Immunoassay
Sample Types	Serum, plasma and other biological fluids.
Uniprot ID	G1T9B3
Sensitivity	0.054 ng/mL
NCBI	XP_002720942
Expiration Date	Please enquire.