


Product Datasheet

Human IL-13/Interleukin-13 ELISA Kit (orb215850)

Description	Human IL-13/Interleukin-13 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G21-N132
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	15816 MW
Uniprot ID	P35225
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1
Expiration Date	Please enquire.



Human IL-13 PicoKine ELISA Kit standard ...