


Product Datasheet

Mouse TIMP-1 ELISA Kit (orb215851)

Description	Mouse TIMP-1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C25-R205
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	22628 MW
Uniprot ID	P12032
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 1024-4380 pg/ml.
Expiration Date	Please enquire.



Mouse TIMP-1 PicoKine ELISA Kit standard...