


Product Datasheet

Mouse MBL2/MPB-C ELISA Kit (orb215854)

Description	Mouse MBL2/MPB-C ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E19-D244
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	25957 MW
Uniprot ID	P41317
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10000, concentration in serum and plasma is 9-60 ug/ml.
Expiration Date	Please enquire.



Mouse MBL2
PicoKine
ELISA Kit
standard c...