


Product Datasheet

Human Kininogen-1 / KNG1 / HMW Kininogen ELISA Kit (orb215856)

Description	Human Kininogen-1 / KNG1 / HMW Kininogen ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q19-S644
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	71957 MW
Uniprot ID	P01042
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is around 600 ng/ml.
Expiration Date	Please enquire.



Human
Kininogen-
1/KNG1
ELISA Kit
standar...