


Product Datasheet

Human IGF2R/IGF-II R ELISA Kit (orb215858)

Description	Human IGF2R/IGF-II R ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S1510-F2108
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	274375 MW
Uniprot ID	P11717
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
IGF2R
PicoKine
ELISA Kit
standard ...