


Product Datasheet

Mouse ANGPTL3 ELISA Kit (orb215859)

Description	Mouse ANGPTL3 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S17-T455
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	52543 MW
Uniprot ID	Q9R182
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 14-39 ng/ml.
Expiration Date	Please enquire.



Mouse
ANGPTL3
PicoKine
ELISA Kit
standar...