


Product Datasheet

Mouse Tissue Factor/F3 ELISA Kit (orb215869)

Description	Mouse Tissue Factor/F3 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A29-E251
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and cell lysates.
MW	32935 MW
Uniprot ID	P20352
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Tissue
Factor/F3
PicoKine
ELISA Ki...