


Product Datasheet

Rat CXCL3 ELISA Kit (orb215873)

Description	Rat CXCL3 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R33-S101
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11109 MW
Uniprot ID	Q10746
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 13-111 pg/ml.
Expiration Date	Please enquire.



Rat CXCL3
PicoKine
ELISA Kit
standard
CU...