

Product Datasheet

HTRA2 ELISA Kit (Human) (orb2158987)

Description	HTRA2 ELISA Kit (Human)
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A134-E458
Storage	Store at 4C for 6 months, at -2°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum and plasma (heparin, EDTA)
Uniprot ID	O43464
Sensitivity	10 pg/ml
NCBI	NP_001308656.1
Expiration Date	Please enquire.