

Product Datasheet

EPB41L1 Antibody (orb215933)

Description

EPB41L1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

E.coli-derived human EPB41L1 recombinant protein (Position: Q691-S881). Human EPB41L1 shares 91% amino acid (aa) sequence identity with both mouse and rat EPB41L1.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for EPB41L1 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

98503 MW

Uniprot ID

[Q9H4G0](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat

Expiration Date

12 months from date of receipt.

100KD —

55KD —

35KD —

25KD —

15KD —

WB analysis of EPB41L1 using anti-EPB41L...



WB analysis of EPB41L1 using anti-EPB41L...