


Product Datasheet

Met (c-Met) Antibody (orb215940)

Description	Met (c-Met) Antibody	
Species/Host	Rabbit	
Reactivity	Human, Rat	
Conjugation	Unconjugated	
Tested Applications	WB	WB analysis of Met (c-Met) using anti-Me...
Immunogen	E.coli-derived human Met recombinant protein (Position: D208-S407). Human Met shares 90% and 91% amino acid (aa) sequences identity with mouse and rat Met, respectively.	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	WB: The detection limit for Met is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	155541 MW	
Uniprot ID	P08581	
Dilution Range	Western blot, 0.1-0.5µg/ml, Rat, Human	
Expiration Date	12 months from date of receipt.	