


Product Datasheet

FOXM1 Antibody (orb215948)

Description	FOXM1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	WB analysis of FOXM1 using anti-FOXM1 an...
Immunogen	E.coli-derived human FOXM1 recombinant protein (Position: A411-Q763). Human FOXM1 shares 75% and 76% amino acid (aa) sequences identity with mouse and rat FOXM1, respectively.	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	WB analysis of FOXM1 using anti-FOXM1 an...
Note	For research use only	
Application notes	WB: The detection limit for FOXM1 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	84283 MW	
Uniprot ID	Q08050	
Dilution Range	Western blot, 0.1-0.5µg/ml, Human	
Expiration Date	12 months from date of receipt.	