



Product Datasheet

ADRA2C antibody (orb216002)

Description	Rabbit polyclonal antibody to Alpha-2C Adrenergic Receptor
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Alpha-2C Adrenergic Receptor. The exact sequence is proprietary.
Target	ADRA2C
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P22086 , P18825 , Q01337
Entrez	24175 , 152 , 11553
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Alpha-2C Adrenergic Receptor.



Immunohistochemical analysis of Alpha-2C Adrenergic Receptor.



Immunofluorescent analysis of Alpha-2C Adrenergic Receptor.