


Product Datasheet

MAP3K8 antibody (orb216065)


Description	Rabbit polyclonal antibody to TPL2 (Phospho-T290)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T290 of human TPL2 protein. The exact sequence is proprietary.
Target	MAP3K8
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q07174 , Q63562 , P41279
Entrez	26410 , 1326 , 116596
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of TPL2 (Phospho-T290)



Immunohistochemical analysis of TPL2 (Phospho-T290)



Immunofluorescent analysis of TPL2 (Phospho-T290)