


## Product Datasheet

### BTBD7 antibody (orb2161)

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| <b>Description</b>         | BTBD7 antibody                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                  |
| <b>Reactivity</b>          | Human, Rat                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                            |
| <b>Tested Applications</b> | FC, IHC-P                                                                               |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human BTBD7 (41-130/1132 aa)              |
| <b>Target</b>              | BTBD7                                                                                   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                   |
| <b>Form/Appearance</b>     | Liquid                                                                                  |
| <b>Concentration</b>       | 1mg/ml                                                                                  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.        |
| <b>Note</b>                | For research use only                                                                   |
| <b>Isotype</b>             | IgG                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                              |
| <b>MW</b>                  | 126 kDa                                                                                 |
| <b>Uniprot ID</b>          | <b>Q9P203</b>                                                                           |
| <b>Dilution Range</b>      | IHC-P=1:100-500, Flow-Cyt=0.2µg/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                         |



Immunohistochemical analysis of formalin...