



Product Datasheet

HNF4A antibody (orb216128)


Description	Rabbit polyclonal antibody to HNF4 alpha (Phospho-S313)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S313 of human HNF4 alpha protein. The exact sequence is proprietary.
Target	HNF4A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P41235 , P22449 , P49698
Entrez	25735 , 15378 , 3172
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HNF4 alpha (Pho...



Immunohistochemical analysis of HNF4 alp...



Immunofluorescent analysis of HNF4 alpha...