



Product Datasheet

PRKCZ antibody (orb216184)


Description	Rabbit polyclonal antibody to PKC zeta
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human PKC zeta. The exact sequence is proprietary.
Target	PRKCZ
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q02956 , P09217 , Q05513
Entrez	5590 , 25522 , 18762
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500, IP: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of PKC zeta expres...



Immunohistochemical analysis of PKC zeta...



Immunofluorescent analysis of PKC zeta
S...