

Product Datasheet

Rabbit IgG (H&L) antibody (orb216231)

Description	Goat polyclonal antibody to Rabbit IgG (H&L)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IH, WB
Immunogen	Rabbit IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000
Expiration Date	12 months from date of receipt.