

Product Datasheet

Human Fab antibody (Biotin) (orb216340)

Description	Mouse monoclonal antibody against Human Fab conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IH, WB
Immunogen	Human IgG Fab
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000
Expiration Date	12 months from date of receipt.