

Product Datasheet

Angiotensinogen peptide (orb216392)

Description	Angiotensinogen peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Arg-Val-Tyr-Ile-His-Pro-Ile-OH
Formula	$C_{43}H_{68}N_{12}O_9$
Purity	> 95%
Source	Synthetic
MW	897.1
CAS Number	102029-49-2
Expiration Date	12 months from date of receipt.