

Product Datasheet

Angiotensinogen peptide (orb216395)

Description	Angiotensinogen peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Sar-Arg-Val-Tyr-Ile-His-Pro-Phe-OH
Formula	$C_{49}H_{71}N_{13}O_{11}$
Purity	> 95%
Source	Synthetic
MW	1018.19
CAS Number	102029-89-0
Expiration Date	12 months from date of receipt.

