

Product Datasheet

Angiotensinogen peptide (orb216402)

Description	Angiotensinogen peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Asp-Arg-Val-Tyr-Val-His-Pro-Phe-OH
Formula	$C_{49}H_{69}N_{13}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1032.18
CAS Number	5649-07-0
Expiration Date	12 months from date of receipt.

