

Product Datasheet

Ac-IETD-pNA peptide (orb216410)

Description	Ac-IETD-pNA peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The chromogenic substrate Ac-IETD-pNA corresponds to one of the cleavage sites of the inactive 32 kD caspase-3 precursor (amino acids 172-175).
Protein Sequence	Ac-Ile-Glu-Thr-Asp-pNA
Formula	$C_{27}H_{38}N_6O_{12}$
Purity	> 95%
Source	Synthetic
MW	638.6
CAS Number	219138-21-3
Expiration Date	12 months from date of receipt.

