

Product Datasheet

Angiotensinogen Human peptide (orb216455)

Description	Angiotensinogen Human peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Arg-Val-Tyr-Ile-His-Pro-Phe-His-Leu-OH
Formula	$C_{58}H_{84}N_{16}O_{11}$
Purity	> 95%
Source	Synthetic
MW	1181.42
CAS Number	56317-01-2
Expiration Date	12 months from date of receipt.

