

Product Datasheet

Atriopeptin III peptide (orb216475)

Description	Atriopeptin III peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Ser-Ser-Cys-Phe-Gly-Gly-Arg-Ile-Asp-Arg-Ile-Gly-Ala-Gln-Ser-Gly-Leu-Gly-Cys-Asn-Ser-Phe-Arg-Tyr-OH (Disulfide bridge Cys7-Cys23)
Formula	$C_{107}H_{165}N_{35}O_{34}S_2$
Purity	> 95%
Source	Synthetic
MW	2549.85
CAS Number	90817-13-3
Expiration Date	12 months from date of receipt.

