

Product Datasheet

Boc-FAAGRK-AMC peptide (orb216498)

Description	Boc-FAAGRK-AMC peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Boc-Phe-Ala-Ala-Gly-Arg-Lys-AMC
Formula	$C_{44}H_{63}N_{11}O_6$
Purity	> 95%
Source	Synthetic
MW	906
Expiration Date	12 months from date of receipt.