

Product Datasheet

Bradykinin (1-3) peptide (orb216506)

Description	Bradykinin (1-3) peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Arg-Pro-Pro-OH
Formula	$C_{16}H_{28}N_6O_4$
Purity	> 95%
Source	Synthetic
MW	368.4
CAS Number	23815-91-0
Expiration Date	12 months from date of receipt.

