

Product Datasheet

Bradykinin Potentiator B peptide (orb216513)

Description	Bradykinin Potentiator B peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	pGlu-Gly-Leu-Pro-Pro-Arg-Pro-Lys-Ile-Pro-Pro-OH
Formula	C ₅₆ H ₉₁ N ₁₅ O ₁₃
Purity	> 95%
Source	Synthetic
MW	1182.46
CAS Number	30892-86-5
Expiration Date	12 months from date of receipt.

