

Product Datasheet

CCAP peptide (orb216528)

Description	CCAP peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Native and synthetic CCAP displayed high ino- and chronotropic effects on a semi-isolated crab heart preparation.
Protein Sequence	H-Pro-Phe-Cys-Asn-Ala-Phe-Thr-Gly-Cys-OH (Disulfide bridge Cys3-Cys9)
Formula	$C_{42}H_{56}N_{10}O_{12}S_2$
Purity	> 95%
Source	Synthetic
MW	957.1
CAS Number	107090-96-0
Expiration Date	12 months from date of receipt.