

Product Datasheet

DAGO peptide (orb216554)

Description	DAGO peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	DAGO, μ -Receptor-specific enkephalin analog.
Protein Sequence	H-Tyr-D-Ala-Gly-N-Me-Phe-Gly-OL
Formula	$C_{26}H_{35}N_5O_6$
Purity	> 95%
Source	Synthetic
MW	513
CAS Number	100929-53-1
Expiration Date	12 months from date of receipt.

