

Product Datasheet

Enterostatin Human peptide (orb216578)

Description	Enterostatin Human peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Ala-Pro-Gly-Pro-Arg-OH
Formula	$C_{21}H_{36}N_8O_6$
Purity	> 95%
Source	Synthetic
MW	496.56
CAS Number	117830-79-2
Expiration Date	12 months from date of receipt.

