

Product Datasheet

Enterostatin, Porcine Rat peptide (orb216579)

Description	Enterostatin, Porcine Rat peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	H-Val-Pro-Asp-Pro-Arg-OH
Formula	C ₂₅ H ₄₂ N ₈ O ₈
Purity	> 95%
Source	Synthetic
MW	582.66
CAS Number	117137-85-6
Expiration Date	12 months from date of receipt.

