

Product Datasheet

Gonadorelin Acetate peptide (orb216627)

Description	Gonadorelin Acetate peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Pyr-His-Trp-Ser-Tyr-Gly-Leu-Arg-Pro-Gly-NH ₂
Formula	C ₅₅ H ₇₅ N ₁₇ O ₁₃
Purity	> 95%
Source	Synthetic
MW	1182.3
CAS Number	71447-49-9
Expiration Date	12 months from date of receipt.

