

Product Datasheet

G-P-R peptide (orb216629)

Description	G-P-R peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	GPR could protect cultured rat hippocampal neurons from Aβ peptide-induced cell death.
Protein Sequence	H-Gly-Pro-Arg-OH
Formula	$C_{13}H_{24}N_6O_4$
Purity	> 95%
Source	Synthetic
MW	328.18
CAS Number	47295-77-2
Expiration Date	12 months from date of receipt.

