

Product Datasheet

Guanylin Human peptide (orb216632)

Description	Guanylin Human peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Endogenous intestinal guanylate cyclase activator.
Protein Sequence	H-Pro-Gly-Thr-Cys-Glu-Ile-Cys-Ala-Tyr-Ala-Ala-Cys-Thr-Gly-Cys-OH (Disulfide bridge Cys4-Cys12, Cys7-Cys15)
Formula	$C_{58}H_{87}N_{15}O_{21}S_4$
Purity	> 95%
Source	Synthetic
MW	1458.68
CAS Number	183200-12-6
Expiration Date	12 months from date of receipt.

